Silicone Super Absorbent Dressing



Healing Phase

1 2 3

Protective Backing

- ▶ Prevents strike through
- ▶ Traps wound exudate in SAP
- ▶ Water proof & bacteria barrier

Super Absorbent Core



Helps minimize maceration risk
Forms gel to build up fluid locking system

Conformable Wicking Layer

- ▶ Rapid exudate wicking -
- ▶ Helps protect fragile tissues
- ▶ Offers patient comfort
- ▶ Will not stick to the wound

OptiSil Silicone Wound Contact Layer -

- ▶ Self-developed silicone oil and perforation technology
- Distributes wound exudate fast and evenly

An excellent dressing for use on heavily exudative wounds. The moisture locking system ensures that exudates are secured within the dressing, reducing the possibility of maceration.

- Excellent exudate management
- Outstanding fluid control
- Non strike-through with bacteria-proof backing
- Minimises the risk of maceration and leakage
- Silicone contact layer for comfort and pain free removal



Please see over for more products







Super Absorbent Dressing



Healing Phase



Absorbency

Protective Backing

- ▶ Prevents strike through
- ▶ Traps wound exudate in SAP
- ▶ Water proof & bacteria barrier

Super Absorbent Core




- ▶ Helps minimize maceration risk-
- Forms gel to build up fluid locking system

Conformable Wicking Layer



- ▶ Helps protect fragile tissues
- ▶ Offers patient comfort
- ▶ Will not stick to the wound





- Excellent exudate management
- Outstanding fluid control
- Non strike-through with bacteria-proof backing
- Minimises the risk of maceration and leakage



Please see over for more products