Vet Pro SP 3000™

Syringe Pump for Special Care

With the Vet Pro SP 3000TM your animals receive the special care they require. Small dose, accurate and linear infusions are readily available as well as automatic brand, size and drug volume detection; making it easy and quick to set up an infusion for any animal. The pump also boasts a rechargeable Li-Polymer battery to allow complete autonomy.

Applications

- Medication/analgesic
- Anaesthetics
- Critical care
- Palliative care
- Chemotherapy
- General fluid replacement

Intended Use

3

- Subcutaneous
- Intravenous
- Percutaneous
- Epidural
- Intra-arterial



Features & Benefits

Automatic Syringe Detection

Automatically detects the syringe brand using 3 independent sensor system.

Unique algorithm detects your drugs volume to reduce dosing errors.

Multi-Program Infusion

Deliver continuous infusion in:

- Lock ON
- Lock OFF
- Rate Mode

Programmable Flow Rates

0.1 - 650 ml/h range to meet your every needs.

Programmable Infusion Time

Instead of rate, you can select from 10 minutes to 100 hours infusion time. The Vet Pro SP 3000™ automatically calculates flow rate for you.

Programmable KVO Rate

0 – 2 ml/h range.

Configure Group Settings

BodyComm[™] software with serial port allows you to customize an infusion protocol and program multiple Vet Pro SP 3000[™] with the same protocol.

Air & Mobile Medicine

CME infusion pumps provide a complete range of infusion onboard any air or ground vehicle. CME pumps are RTCA-D0160 compliant.

BodyComm[™] Software

Review patient history with charts and 500 events log or manage pump settings on your PC.

True Portability

Size & Weight

205 x 98 x 45 mm 411 g

Battery-Powered

Uses a long-life rechargeable Li-Polymer battery for up to 100 full syringe deliveries.

Flexible Configuration

Human Pump Simplified

Software interface optimized for veterinary use.

Supports Multiple Languages

Over 20 languages available to meet your needs and preferences.

Multi Pump Connection

MultiDock[™] also enables connection and charging of several pumps together with single power cable.

PCA Model Available

Designed especially to deliver analgesics using a Luer-Lock syringe.

Battery or AC-Powered

Change between battery and AC (mains) power without interrupting infusion therapy.

Safety & Security

Dose Titration

Operator or patient can titrate flow rate up to pre-set limit for safe, convenient dose adjustment.

Full Alarm Protection

Helps prevent events and instill confidence in veterinarians.

Easy-Scroll Event Review

Press INFO key for hour-by-hour event log.

3-Level Access Codes

Prevent accidental or unauthorized program changes.