

Hydrogel Dressing



Healing Phase

1

2

3

4

Hydration



Hydrogel


- ▶ Hydrophilic polymers
- ▶ Donates water
- ▶ Comfortable and flexible

Blue Dots Mesh

- ▶ Backbone of the dressing
- ▶ Easy for wound assesment

Protective Backing

- ▶ High MVTR (Moisture Vapour Transmission Rate)
- ▶ Smooth and conformable
- ▶ Waterproof & bacteria barrier



For minimal exudative and dry wounds. The dressing has the ability to absorb a minimum amount of fluid by swelling, but also has the ability to donate moisture to the wound. The hydrogel dressing is constructed of a semi-transparent flexible gel. This cool dressing provides optimal conditions for tissue granulation and epithelisation.

- Provides a continual moist environment for wound healing
- Rehydrates and loosens necrotic tissue
- Optimal conditions for granulation and epithelisation



Please see over for more products

01347 878697

www.jakmarketing.co.uk

Amorphous Hydrogel



Healing Phase

1

2

3

4

Amorphous Hydrogel

- ▶ Gently re-hydrates necrotic tissue
- ▶ Facilitates autolytic debridement
- ▶ Creates a moist wound environment



Hydration



A colourless, cooling, transparent gel with excellent performance. The non-adherent gel donates moisture to the wound and maintains an optimal healing environment to facilitate debridement and epithelisation.

- Cooling and soothing on application
- Easy to apply to all wounds
- Crystal clear transparency for easy visualisation of wound assessment
- Maintains hydration balance to provide an optimal healing environment
- Encourages autolytic debridement



Please see over for more products

01347 878697

www.jakmarketing.co.uk